



Sharps Passing Tray

Secure and safe transfer of Sharps

The innovative Sharps Passing Tray reduces the potential risk of needlestick injury to clinicians.

Globally more than 2 million¹ occupational exposures to sharp injuries occur among 35 million healthcare workers each year. Specifically, 18,000² in Australia, 385,000³ in the US, and 1,000,000⁴ in the EU.

The versatile design improves safety with ergonomically designed textured grip on sides, top and base while a channel at each end secures variable sharps component sizes up to a 30mL Syringe facilitating safer transfer.

The low profile design maximises limited storage space and reduces biohazard waste.

Partners in Healthcare Associated Blood Stream Infection Prevention Products and Strategies



¹https://www.who.int/occupational_health/topics/needleinjuries/en/

²<https://www.mtaa.org.au/sites/default/files/uploaded-content/website-content/Sharpsv5.pdf>

³<https://www.cdc.gov/sharpsafety/>

⁴<https://www.europeanbiosafetynetwork.eu/wp-content/uploads/2017/01/EU-Sharps-Injuries-Implementation-Guidance.pdf>

Sharps Passing Tray

- Healthcare workers have a duty to use a hands free or neutral zone, i.e. an agreed upon location where sharps are placed for transfer so no two individuals are touching the sharp simultaneously¹.
- ITL's Sharps Passing Tray patented versatile design enables a safer needle and sharps transfer between clinicians by shielding the sharp section of the instrument while leaving the handle free to access.
- Sharps Passing Tray accommodates all commonly used sharps used in the operating room making healthcare workplaces safer
- Clinically disposable and made from lightweight, durable medical grade polypropylene suitable for recycling

¹ ACORN Statement: Management of Sharps in the Perioperative environment 1.3.2



Small

Large

Protects

- Three different pick-up zones enable the careful placement and removal of components
- An innovative channel at each end secures the syringe flange for safe transfer

Ergonomic & Stable

- Wide base provides optimal stability on uneven surfaces
- Tactile grip surfaces on both sides, top and base helps to negate any handling slippage

Maximise Efficiencies

- Low profile optimizes storage
- Available in two sizes and can be easily added to customised procedure packs to reduce overall packaging

For further information or samples, visit itlbiomedical.com or contact your local sales office.

Product Codes

A500103	Sharps Passing Tray Small Sterile	104 / Carton
A500102	Sharps Passing Tray Small Non Sterile	115 / Carton
A500117	Sharps Passing Tray Large Sterile	50 / Carton
A500116	Sharps Passing Tray Large Non Sterile	75 / Carton

Product Information



www.itlbiomedical.com/symbolglossary