



THERE'S ONLY ONE SAFETY CAP THAT
SECURELY FITS BOTTLES AND TEST TUBES.

SampLok[®] Adapter Cap

Safer Blood Culture Sample Collection

INTEGRATED SAFETY for collecting blood
samples for blood culture testing

REDUCES RISK of inadvertent
needlestick injury

FLEXIBLE DESIGN allows sampling with
a variety of culture bottles and vacuum
test tubes

www.itlbiomedical.com



SampLok® Adapter Cap

Safer Blood Culture Sample Collection

Samplok® Adapter Cap is a versatile device for safe blood sample collection into blood culture bottles and test tubes.

Samplok® Adapter Cap fits securely over the blood culture bottle top.

Fits a variety of culture bottles and vacuum tubes.

PRODUCT FEATURES

The cap ensures the sampling needle is well positioned in the culture bottle septum.

At the end of the blood sample collection the lid closes to shield the needle away.

Available with or without a pre-attached luer.



Optional insert for use with sample tubes and long neck culture bottles.



Can be attached to blood collection and winged needle sets. The integrated safety lid provides added protection from the sampling needle.



Product Information

 Latex Free

 Single Use Only



Product Codes

A100706 Without luer, non-sterile, 350/case

A100710 With luer, individually packaged, sterile, 100/case

A100709 Insert, non-sterile, 600/case, 50/box, 12 boxes/case



HUMAN INSPIRED. SCIENTIFICALLY ENGINEERED.

ITL BioMedical

www.itlbiomedical.com

Sales Office: 1.888.411.2851

sales@itlbiomedical.com