

A SEGMENT SAMPLER
OFFERING SAFETY & MEASURED SAMPLES
ALL IN ONE UNIQUE DEVICE

SampLok[®] Segment Sampler

SIMPLE – Efficient & safe sampling procedure

UNIQUE PROCESS – Obtain a measured sample from segment tubing

ERGONOMIC ADVANTAGES – Reduced biohazard exposure and cross contamination risk

www.itlbiomedical.com



SampLok[®] Segment Sampler

MEASURED SAMPLE SAFE AND FAST

Now there is a safer way to collect measured samples from tubing segments including

- Platelet samples from platelet bag tubing segments
- Component samples from storage bag integral tubing segments



PRODUCT FEATURES

Recessed sampling needle in the Segment Sampler provides a safe way to collect blood component samples from tubing segments

Reduces the risk of sharps injury

Measured samples with needleless syringe



Attaching a needleless syringe with volume graduations to the Segment Sampler makes this the safe segment sampling device that provides a measured sample.

ADDITIONAL ADVANTAGES

Ergonomic process enhancements and reduced biohazard exposure – Segment Sampler eliminates cutting and squeezing of tubing segments as well as dripping sample from an open segment.



Product Information

-  Latex Free
-  DEHP Free
-  Single Use Only



Product Code

A200200 Individually packaged, sterile, 360/case



HUMAN INSPIRED. SCIENTIFICALLY ENGINEERED.

ITL BioMedical

www.itlbiomedical.com
Sales Office: 1.888.411.2851
sales@itlbiomedical.com