



## Safety SubCulture Unit 2

### Versatile dispensing unit for a safer subculture technique

Make your lab safer and more efficient with Safety SubCulture Unit 2.

It enables simple and safe needleless transfer of samples from culture bottles to slides and media plates.

Designed to fit securely onto a range of culture bottles, the plastic dispensing tip makes it easy to dispense samples onto slides and plates.

The protective cap captures gas and culture discharge, reducing potential exposure.

The new standard for safety in blood culture sampling



# Choose the safer solution for subculture

## Safety SubCulture Unit 2:

- Needleless dispensing of blood culture samples, reducing the risk of sharps injury
- Secure fit onto culture bottle septum
- Protective cap over the dispensing tip reduces potential exposure
- Safe and easy sample dispensing procedure onto slides and plates
- Fits a wide variety of culture bottles



Visit [itlbiomedical.com](http://itlbiomedical.com) to request your sample today

### Product Code:

A100722: Sterile, 288/case, 48/box, 6 boxes/case

### Product Information:



Not Made with Natural Rubber Latex



Single Use Only

